




Gradient Segmentation

Segmentation

- To quantify cell shape and motion, cells must be identified in 3D and then tracked over time
 - Velocity
 - Distance
 - Dendrite Extension
 - Cell-Cell interactions
- Many programs exist to do this, with varying levels of success




Imaris Algorithm

- On 1 channel, set a threshold separating foreground and background
- Smooth off corners, filter small objects
- Watershed transform to separate touching objects
 - Relies on indentation points




Problems with Imaris Segmentation

- Assumes spherical cells that are slightly touching
- Assumes sharp change in intensity on cell boundaries









- Aortic Macrophages:
 - Low resolution
 - 0-1 pixels between cells
- Many true double positives, with small difference in ratio of fluorophores
 - Thresholding on GFP or YFP will always pick up both cells
 - Using GFP/YFP or YFP/GFP ratio is very noisy and condenses all double positives into a small range
- Dendritic-shaped with indentation points that do not correlate to cell boundaries





Imaris Segmentation- Giant Errors



Errors



Errors



Macrophage Segmentation

- Often the boundary between cells isn't defined by a sharp change in GFP or YFP intensity, but a change in the apparent color
- Need a segmentation algorithm that:
 - Uses both GFP and YFP information simultaneously (NOT the ratio)
 - Blind to shape

Dynamic Color Gradient Segmentation

Unsupervised color image segmentation using a dynamic color gradient thresholding algorithm

Guru Prashanth Balasubramanian^{*a}, Eli Saber^a, Vladimir Misic^b, Eric Peskin^a, Mark Shaw^c

Human Vision and Electronic Imaging XIII, edited by Bernice E. Rogowitz, Thrasyvoulos N. Pappas,
Proc. of SPIE-IS&T Electronic Imaging, SPIE Vol. 6806, 68061H, © 2008 SPIE-IS&T · 0277-786X/08/\$18

Original



Segmented



Dynamic Color Gradient Segmentation

- Uses multiple color gradients simultaneously
- Includes “texture” (measure of variance within a cell)
 - Dim cells have more “texture” than bright cells
- My adaptations:
 - 3D instead of 2D
 - Does not process background pixels to greatly increase efficiency
 - Multiple rounds of processing to better incorporate known facts about cells

Algorithm Overview








Fig. 2. Flowchart of the proposed algorithm.

1- Remove Background

- In Imaris, make surfaces over GFP+YFP channel to cover cell volume
- Mask those pixels and force all background pixels to be 0
- Import the cleaned GFP and YFP data into Matlab




2- Calculate 2D dual color gradient maps

- For each pixel calculate the total gradient in both colors in 2D
- Set all background pixels to 0 without calculating



3- Enhance the Gradient Map

- Using a combination of high, low, and middle thresholds, increase the contrast on the gradient map
 - Some are set automatically, some are set manually and will need to be tweaked




4- Seed generation

- Find regions below a threshold on the 2D gradient map (even color in raw data)
- These are the starting seeds of region growing
 - Size and threshold can be adjusted



5- Region Growing

- Increase the max threshold on the gradient map
- For each newly selected area, decide if it should be merged with an existing region or start it's own region
 - Based on size and comparison of median intensity of GFP and YFP



6- Region Merging Round 1

- Region growing results in a highly oversegmented image (~3000 regions / timepoint in this example)
 - If neighboring cells are merged here, stop and change some variables
- **In 3D**, make a graph of all adjacent neighbor regions
- For each touching pair, compute a similarity measure:


$$S_{A,B} = (\mu_A - \mu_B)^T * inv(\text{COV}_A + \text{COV}_B) * (\mu_A - \mu_B)$$

Difference in mean
intensities in both channels

Texture, or covariance
between GFP and YFP

Lower means more similar


- Result is a list of similarities measurements between touching regions



6- Region Merging


- Merge the most similar regions (balance efficiency vs accuracy)
- Iteratively repeat calculating connections graph and similarity measurements until smallest similarity measure is above a threshold
 - Most sensitive variable in the program

- Regions are now in 3D




7- Smart Region Clean Up

- Forcibly merge regions to most similar touching neighbor that fail certain criteria:
 - Minimum size
 - Landlocked- Completely enclosed by one or more cells
 - Does not apply to all potential imaging, but does work for atherosclerotic macrophages





8- Region Merging Round 2

- Repeat merging procedure on “cleaned up” regions to the same threshold




9- Form Surfaces

- For each identified region, import the intensity of the voxels back into Imaris, and have Imaris create a surface around it
 - Smoother than having Matlab do it




Improvements with this algorithm


Raw




Imaris Default




Color Gradient




Raw



Imaris Default




Color Gradient




Improvements


Raw




Imaris Default




Color Gradient




Raw



Imaris Default




Color Gradient




Errors in this algorithm

Raw



Imaris Default



Color Gradient

